

Instructions for use of UT block

Chemical composition

Components	UT
ZrO ₂ +HfO ₂ +Y ₂ O ₃	≥99%
Y ₂ O ₃	9%-10%
Al ₂ O ₃	<0.05
Other Oxides	≤0.05

Physical characteristics

Properties	UT
Density before sintering (g.cm ⁻³)	3.15-3.25
Density after sintering (g.cm ⁻³)	6.05-6.07
CTE(25-500°C)(K ⁻¹)	10.5
Flexural strength after sintering(Mpa)	650
Accelerated aging surface monoclinic phase content	<15%
Light transmittance	49%
Chemical solubility after sintering(μ g.cm ⁻²)	<100
Cytotoxicity	0 Level
Radioactivity(Bq.g ⁻¹)	<0.1
Sintering temperature(°C)	1450-1530

Instructions for use

Bloomden UT (ultra translucent) can be used for full contour and coping crown. The following instructions provide general guidelines for designing, milling, coloring, sintering, polishing and glazing and should be followed very carefully to avoid any loss of aesthetics, fit, durability or overall quality.

Designing Bloomden UT

Design Option	Design Guidance
Drill Compensation	Drill compensation must be activated for all substructures milled from a solid structure.
Cement Gap	The distance where the coping intersects the die at the margin area. Use this setting to control margin fit.
Extra cement Gap	This distance between the coping walls and the die. Use this setting to control internal fit.
Distance to Margin Line	The distance from the margin outer line to the start of the interior

	wall of the coping.
Smooth Distance	The distance from the margin line to the margin engagement point should be set at 0.2mm.
Drill Radius	The drill radius is the size of the smallest end mill used to mill the pattern.
Drill Compensation Offset	The distance from the margin line to the area affected by drill compensation should be a minimum of 0.5mm.
Margin line Offset	The effective thickness of the margin line and should not be less than 0.16mm. Thinner margin lines will result in a higher failure rate.
Offset Angle#1	The offset angle should not be less than 65°.
Extension Offset	The extension offset should be less than 0.01mm.
Wall thickness	A nominal wall thickness of 0.5mm will ensure a consistently quality product. Reducing this value could result in fractures or holes in the framework.
Bridge Connectors	Recommended Anterior restorations: 6mm ² minimum. Recommended Posterior restorations: 9mm ² minimum.

Milling Bloomden UT

Pre-sintered zirconia material has an inherent shrinkage rate associated with each production. This number of shrinkage rate must be input into the milling preparation software to ensure the accuracy of the eventual restoration.

When milling Bloomden UT, always follow these general guidelines:

- Only use sharp end mills with carbide or diamond coating.
- Do not use any restoration that has chips or cracks. Remove the units from the disc using a handpiece with a diamond-coated burr.
- Smooth the support areas with a medium-grit rubber polishing wheel.
- Remove any residual zirconia dust with an art brush.
- If a wet mill is used make sure all the zirconia is completely dry before sintering. Air dry for at least 15 minutes prior to sintering. Damp zirconia will crack if placed in the sintering oven.

Coloring Bloomden UT

Bloomden UT is compatible with Bloomden UT- A1, A2, A3, B1, B2, B3 6 shades coloring liquid.

Sintering Bloomden UT

Temperature rising	Temperature	Time span	Total time
<u>°C/min</u>	<u>°C</u>	<u>hour</u>	<u>hour</u>
0	20		0
10	900	1.5	1.5
0	900	0.5	2
3.27	1450	2.8	4.8
0	1450	2	6.8
4.58	900	2	8.6

Sintering temperature of Bloomden UT is from 1450°C to 1530°C, you should make sure peak temperature at this range.

Polishing and glazing Bloomden UT

Use Bloomden polishing tool specially made for zirconia to adjust the occlusion and shape of the restoration. After coarsely finishing, proceed to external dyeing and glazing.

Finished

The production of restoration is completed.